



HotDog Veterinary Warming System

HotDog warming represents a significant breakthrough in technology and features.

The patented ThermAssure™ fabric's conductive polymer responds to low-voltage, direct current by becoming uniformly warm without hot spots.



Benefits :

- Quiet.
- Fast Warm-up time.
- Small and lightweight.
- Even-Heat Distribution.
- Environmental friendly.
- 7-inch touch screen operation.
- Varies size, water proof blankets.
- Great for surgery, dental procedure, x-ray, and recovery.



For more information, Please contact Asia Vet Medical Limited.



Asia Vet Medical Limited
Tel: (852) 2335 9200
Fax: (852) 2335 0008
Email: sales@avmhc.com